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
APPLICATION GUIDE

Application Studies (grouped according to product type) are presented as a tool for beginning research in relevant product sectors. They assist in providing basic methods, explanations and tips for commencing new tests or investigating new products and will hence be subject to change according to specific requirements and product differences.

MATERIALS & PRODUCTS	FOOD
ADHESIVES	BAKERY
COSMETICS & SKINCARE	CEREALS
PHARMACEUTICAL & MEDICAL	CONFECTIONERY
PHARMACEUTICAL PACKAGING	DAIRY
PERSONAL CARE PRODUCTS	FRUIT
MEDICAL DEVICES	VEGETABLES
CONTROLLED RELEASE	FISH
HAIR & HAIR PRODUCTS	MEAT
ELECTRONICS	PASTA & RICE
LEATHER	SNACKS
MATERIALS	FOOD PACKAGING
GELS & FILMS	PETFOOD & ANIMAL FEED
POWDER & GRANULES	GELS & FILMS

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TEST METHODS FOR GELS & FILMS

Click here for an overview of Testing Attachments for Gels & Films

See what others have published in this Application Area

A SELECTION OF APPLICATION STUDIES - TEST METHODS & ANALYSIS

Air Freshener Gel
☰ [Determination of gel strength by penetration with a cylinder probe \(air freshener gel\)](#)

Breath Strips
☰ [Comparison of the film burst strength of fresh and 13 day old breath strips](#)

Drug Release Films
☰ [Comparison of the film burst strength of two pharmaceutical drug release films](#)
☰ [Comparison of the film relaxation properties of two pharmaceutical drug release films](#)
☰ [Comparison of the film resilience properties of two pharmaceutical drug release films](#)

Gels
☰ [Comparison of rupture force and brittleness/elasticity of three gel formulations as a means of monitoring effects of quality, concentration and processing](#)

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Product:	GEL
Objective :	Comparison of rupture force and brittleness/elasticity of three gel formulations as a means of monitoring effects of quality, concentration and processing

Investigation Purpose:


The measurement of gel strength is of widespread interest in the manufacture of pharmaceutical, medical and cosmetic product: pectin, gelatin, agar etc. are important in the development of such products as coronary stents where hydrogel polymers are self strong, superficial resemblance to living, soft tissue. Other products for which gel-forming properties are useful are in the manufacture of suppositories, soft gel capsules and bacterial growth media. The strength of gels can also be utilised in products such as toothpaste required end product.

Standard probes such as that required for gelatine testing (BS 757), as used in this procedure, are available for the assessment of BI a cylinder probe which possesses a slightly rounded (or termed 'radiused') edge which makes the probe less sharp and therefore d the probe, which is often the case for standard cylinder probes when testing gels.


TA Settings :	Mode:	Measure Force in Compression
	Option:	Return to Start
	Pre-Test Speed:	1.0 mm/s
	Test Speed:	0.5 mm/s
	Post-Test Speed:	10.0 mm/s
	Distance:	8mm
	Trigger Type:	Auto - 10g
	Tare Mode:	Auto
	Data Acquisition Rate:	400pps
	<input type="button" value="LOAD PROJECT"/>	

Accessory: [0.5" Radius Cylinder \(P/0.5R\)](#) using 5kg load cell

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