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The screenshot displays the 'APPLICATION GUIDE' window. On the left is a 'CONTENTS' sidebar with a tree view of various product categories. The main area shows a grid of two columns: 'MATERIALS & PRODUCTS' and 'FOOD'. Each category in the grid is represented by a small image and a text label. A red arrow points to the 'PASTA & RICE' category in the 'FOOD' column. Below the grid is a section for 'INTERNATIONAL STANDARD METHODS'. The right side of the window shows a blank workspace with a 'None' label and a 'Time (sec)' scale at the bottom.

MATERIALS & PRODUCTS	FOOD
ADHESIVES	BAKERY
COSMETICS & SKINCARE	CEREALS
PHARMACEUTICAL & MEDICAL	CONFECTIONERY
PHARMACEUTICAL PACKAGING	DAIRY
PERSONAL CARE PRODUCTS	FRUIT
MEDICAL DEVICES	VEGETABLES
CONTROLLED RELEASE	FISH
HAIR & HAIR PRODUCTS	MEAT
ELECTRONICS	PASTA & RICE
LEATHER	SNACKS
MATERIALS	FOOD PACKAGING
GELS & FILMS	PETFOOD & ANIMAL FEED
POWDER & GRANULES	GELS & FILMS

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Click here for an overview of Testing Attachments for Pasta Products

See what others have published in this Application Area

A SELECTION OF APPLICATION STUDIES - TEST METHODS & ANALYSIS

Gnocchi

- [Firmness comparison of two brands of gnocchi by compression](#)

Lasagne

- [Comparison of breaking stress/strength of 3 types of dry lasagne using a three-point bend test](#)
- [Comparison of pasta stickiness cooked in hard and soft water](#)
- [Assessment of the tensile strength and extensibility of three types of cooked lasagne](#)


Noodles

- [Comparison of elasticity \(or 'tensile strength'\) of noodles](#)
- [Comparison of hardness and adhesiveness of noodles using a cylinder probe](#)


8. คลิ๊ก LOAD PROJECT

Product:	NOODLES	
Objective :	Comparison of elasticity (or 'tensile strength') of noodles	
TA Settings :	<i>Mode:</i>	Measure Force in Tension
	<i>Option:</i>	Return To Start
	<i>Pre-Test Speed:</i>	1.0 mm/s
	<i>Test Speed:</i>	3.0 mm/s
	<i>Post-Test Speed:</i>	10.0 mm/s
	<i>Distance:</i>	100mm
	<i>Trigger Type:</i>	Auto - 5g
	<i>Tare Mode:</i>	Auto
	<i>Data Acquisition Rate:</i>	200pps
	LOAD PROJECT	
Accessory:	Spaghetti tensile grips (A/SPR) using 5kg load cell	

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