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
Contents Index Search

CONTENTS

- Test Methods by Product T
- Test Methods According to
- Guide to Measuring a Spec
- Introduction to Texture Ana
- Data Analysis Techniques
- Advanced Techniques
- Probes & Fixtures Guide
- Measure More Than Force
- Sample Macros
- Sample Projects
- Sample Sax Basic Scripts
- Published papers using the
- Video Clip Library
- Tutorials & User Tips
- Frequently Asked Question
- Library Tests Available
- View YouTube Channel
- Open Website
- Request Help
- Request a Brochure
- Request a Manual or Post

APPLICATION GUIDE

Application Studies (grouped according to product type) are presented as a tool for beginning research in relevant product sectors. They assist in providing basic methods, explanations and tips for commencing new tests or investigating new products and will hence be subject to change according to specific requirements and product differences.

MATERIALS & PRODUCTS	FOOD
ADHESIVES	BAKERY
COSMETICS & SKINCARE	CEREALS
PHARMACEUTICAL & MEDICAL	CONFECTIONERY
PHARMACEUTICAL PACKAGING	DAIRY
PERSONAL CARE PRODUCTS	FRUIT
MEDICAL DEVICES	VEGETABLES
CONTROLLED RELEASE	FISH
HAIR & HAIR PRODUCTS	MEAT
ELECTRONICS	PASTA & RICE
LEATHER	SNACKS
MATERIALS	FOOD PACKAGING
GELS & FILMS	PETFOOD & ANIMAL FEED
POWDER & GRANULES	GELS & FILMS

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Click here for an overview of Testing Attachments for Snack Products

See what others have published in this Application Area

A SELECTION OF APPLICATION STUDIES - TEST METHODS & ANALYSIS

Aerated Desserts

- [To measure the acoustic emissions and textural properties during back extrusion of aerated desserts](#)


Crisps & Chips

- [Comparison of textural qualities of crisp samples by bulk compression using an Ottawa cell](#)
- [Fracturability of Tortilla Chips](#)
- [Measurement of fracturability of an extruded snack sample using a 5-bladed Kramer shear cell](#)
- [Firmness measurement of chips/french fries with multiple chip rig](#)


8. คลิ๊ก LOAD PROJECT

Product:	TORTILLA CHIPS	
Objective:	Fracturability of Tortilla Chips	
TA Settings :	<i>Option:</i>	Measure Force in Compression
	<i>Pre-Test Speed:</i>	1.0 mm/s
	<i>Test Speed:</i>	1.0 mm/s
	<i>Post-Test Speed:</i>	10.0 mm/s
	<i>Distance:</i>	3mm
	<i>Trigger Type:</i>	Auto - 5g
	<i>Tare Mode:</i>	Auto
	<i>Post-Test Wait Time:</i>	500pps
	LOAD PROJECT	
Accessory:	Crisp Fracture Rig (HDP/CFS) using 25kg load cell Heavy Duty Platform (HDP/90)	

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