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
APPLICATION GUIDE

Application Studies (grouped according to product type) are presented as a tool for beginning research in relevant product sectors. They assist in providing basic methods, explanations and tips for commencing new tests or investigating new products and will hence be subject to change according to specific requirements and product differences.

MATERIALS & PRODUCTS	FOOD
ADHESIVES	BAKERY
COSMETICS & SKINCARE	CEREALS
PHARMACEUTICAL & MEDICAL	CONFECTIONERY
PHARMACEUTICAL PACKAGING	DAIRY
PERSONAL CARE PRODUCTS	FRUIT
MEDICAL DEVICES	VEGETABLES
CONTROLLED RELEASE	FISH
HAIR & HAIR PRODUCTS	MEAT
ELECTRONICS	PASTA & RICE
LEATHER	SNACKS
MATERIALS	FOOD PACKAGING
GELS & FILMS	PETFOOD & ANIMAL FEED
POWDER & GRANULES	GELS & FILMS

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Packaging Seals


- Measurement of seal strength of ribbed heat-sealed foil packaging
- Measurement of peel strength of a container lid seal using a Peel Strength Test Rig
- Assessment of the peeling characteristics of packaged sliced meat using Articulated Tensile Grips

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
Product:	PACKAGING SEAL	
Objective:	Measurement of seal strength of ribbed heat-sealed foil packaging	
TA Settings:	<i>Mode:</i>	Measure Force in Tension
	<i>Option:</i>	Return To Start
	<i>Pre-Test Speed:</i>	1.0 mm/s
	<i>Test Speed:</i>	1.5 mm/s
	<i>Post-Test Speed:</i>	10.0 mm/s
	<i>Distance:</i>	32mm (at least twice the seal length)
	<i>Trigger Type:</i>	Auto - 5g
	<i>Tare Mode:</i>	Off
	<i>Data Acquisition Rate:</i>	400pps
	LOAD PROJECT	

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The screenshot shows the Exponent software interface with a data table. The table has columns A through F. Column A contains test IDs and descriptions. Column B contains batch names. Column C contains formulas. Column D contains Rib Seal Strength values. Column E contains Force 1 values. Column F contains Distance 1 values. The table is divided into sections for 'walkerseal' and 'film food r' tests, each with an 'Average', 'S.D.', and 'C.V.' row. The 'End of Test Data' row is highlighted in purple.

	A	B	C	D	E	F
1	Test ID	Batch		Rib Seal Strength N/mm	Force 1 kg	Distance 1 mm
2				Strength 1	Force 1	Distance 1
3						
5	Start Batch walkerseal	walkerseal				
6	Packaging Seal Strength	walkerseal		.091		
7	End Batch walkerseal	walkerseal				
8	Average	walkerseal (F)	AVERAGE("BATCH")	.091		
9	S.D.	walkerseal (F)	STDEVP("BATCH")	.000		
10	C.V.	walkerseal (F)	STDEVP("BATCH")/AVERAGE("BATCH")*100	.000		
11	Start Batch film food r	film food r				
12	film food r5	film food r			.347	40.325
13	film food r6	film food r			.236	23.555
14	film food r7	film food r			.267	34.235
15	End Batch film food r	film food r				
16	Average	film food r (F)	AVERAGE("BATCH")		.283	32.705
17	S.D.	film food r (F)	STDEVP("BATCH")		.047	6.931
18	C.V.	film food r (F)	STDEVP("BATCH")/AVERAGE("BATCH")*100		16.549	21.193
19	End of Test Data					

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